



HOW TO RECOGNIZE AN ORIGINAL EVOLIS HIGH TRUST® CONSUMABLE



Hereafter is a description of the several points to check depending on your ribbons.

PRIMACY, PRIMACY2 & ZENIUS RIBBONS IN CASSETTE (P/N: RXXX)



- Check the **Evolis High Trust®** and **Evolis logos** on the label.



- Two 'O' from the Evolis logo are placed on each plastic ribbon core.



Cores made of POM material



Cores made of rPET material


PRIMACY & PRIMACY2 RIBBONS WITHOUT CASSETTE (EASY4PRO RANGE, P/N: NXXX)

- We placed two 'O' from the Evolis logo on each plastic ribbon core.



RETRANSFER RIBBONS: AVANSIA/AGILIA (P/N: RTXXX)

- Check the **Evolis High Trust®** and **Evolis logos** on the labels.



- Check the **Evolis logo** on the RFID tag for Avansia ribbons



COLOR RIBBONS WITH PART NUMBERS R3012/R3011/R3511/R3013/R3514 IN INDIVIDUAL BOXES (FOR SECURION AND QUANTUM2 PRINTER):

- Check the Evolis High Trust® and Evolis logos on the box.
- We placed a **holographic sticker** on the ribbon box, with several security features, such as a rainbow color fine lines guilloche design, nano text, two channel images, etc.
- We also integrated a unique feature: **the floating image**. To check the authenticity of your ribbon, you just have to pass the sticker under a standard light and the 'O' from the Evolis logo will appear.



- On each side of the ribbon cores, check the logo



MONOCHROME RIBBONS WITH PART NUMBERS R2011/R2131 IN SEALED PLASTIC BAGS (FOR SECURION AND QUANTUM2 PRINTERS):

- Check the Evolis High Trust® and the Evolis logos on the label.
- On each side of the ribbon cores, check the logo



Evolis - 14 avenue de la Fontaine - Z.I. Angers-Beaucouzé 49070 Beaucouzé - France
T +33 (0) 241 367 606 - F +33 (0) 241 367 612 - info@evolIS.com

www.evolIS.com